

GelView Master

Gel Documetation System

Gel in Touch



IPad or Android tablet



PC



Wi-Fi connection for wire-free lab environment



- Easy to use with quick imaging capture
- Intuitive imager software with USB connection to PC
- Multi-image file saving selection: BMP, TIFF, JPEG and PNG
- Built-in, 21 x 26 cm pullout UV transilluminator for ease of use and gel cutting
- Blue light technology for operation and environment safe lab experiment
- Pre-focused camera lens, no adjustment needed
- Optional blue LED or white light table available, easy installation
- Ready to use right out of the box
- Imager analysis software included

GelView Master
Gel Documetation System





GelView Master (standard version)



GelView Master Wi-Fi (Wi-Fi version)

GelView Master Imager Systemt Specification

Model	GelView Master (standard version)	GelView Master Wi-Fi (Wi-Fi version)
Camera / Lens	5.0 MP scientific grade camera	5.0 MP scientific grade camera
Max. Aperture	f/1.2	f/1.2
Built-in UV transilluminator	Yes, 312 nm (optional white or blue light available)	Yes, 312 nm (optional white or blue light available)
Image Storage	PC only	PC / iPad / Android tablet
Safety Device	Safety door switch	Safety door switch
Features	USB connects to a PC Image Capture software included Connect to PC; PC required	Wifi to PC / iPad / Android tablet Image Capture software included Connect to PC / iPad / Android tablet PC / iPad / Android tablet required

GelView Master Ordering Information

PRODUCT	CATALOG NUMBER#
GelView Master Imaging System, standard version	GD-GVMS
GelView Master Imaging System, standard version, UV high performance	GD-GVMS-UV
GelView Master Wi-Fi Imaging System, wifi version with iPad compatibility	GD-GVMW
GelView Master Wi-Fi Imaging System, wifi version with iPad compatibility, UV high performance	GD-GVMW-UV

Headquarters

Dynamica Scientifi Ltd.
 Unit 12 Cromwell Business Centre,
 Howard Way, Newport Pagnell,
 MK16 9QS, United Kingdom
 P: +44(0)1908 211 900
 F: +44(0)1908 211 909
 Email: info@dynamica-eu.com

Australasia

Dynamica Pty Ltd
 Unit C3 Hallmarc Business Park
 Corner Centre and Westall Roads
 Clayton, Victoria, 3168, Australia
 P: +61 3 8540 5988
 F: +61 3 9548 7177
 Email: info@dynamica.com.au

Asia

Dynamica (Asia) Limited
 6/F., Mita Ctr.,552-566 Castle
 Peak Rd.,Kwai Chung, Hong Kong
 P: +852 3583 1581
 F: +852 3583 1580
 Email: info@dynamica-asia.com
 Web: www.dynamica-eu.com