



SERIES 390

Release software version 1.0.5 and launch of the density kit

HIGHLIGHTS OF THE RELEASE 1.0.5

- Density determination application
free for HF and HA models
optional for HE and HM models (free evaluation); 390-8609
 - guided workflow with animated step by step instructions for simple and user friendly operation
 - supports solid and porous materials as well as liquids
 - supports density determination with the density kit
 - supports density determination „below balance“
 - fully integrated context help
 - integrated installation instructions for the „density kit“ and „below balance“
 - allows reliable and error-free work
- After calibration a report is automatically generated
- After adjustment and levelling return directly to the previous app via home-button



LAUNCH OF THE DENSITY KIT SERIES 390

For easy handling in density determination, we recommend the density kit. It is fully supported by the density application.

ITEM NUMBERS

350-8719: Density kit

350-8720: Density kit solids only (without glass body and hook)

* Glass body for density determination of liquids

